permobil | SEATING + POSITIONING

ROHO® Hybrid Elite Cushion

Operation Manual

s h a p e f i t t i n g t e c h n o l o g y®

Contact Information

 ROHO, Inc.
 100 N. Florida Avenue
 Belleville, IL 62221-5429 U.S.A.
 U.S.A.:
 800-851-3449
 Fax 888-551-3449
 Customer Support:
 orders.roho@permobil.com
 Outside the U.S.A.:
 618-277-9150
 Fax 618-277-6518
 intlorders.roho@permobil.com
 permobilus.com

 **Warning**

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in serious injury.

 **Caution**

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in minor to moderate injury or damage to the equipment or other property.

Incident Reporting

In the event of a serious incident or injury, contact Customer Support.

Contents

Intended Use			2
Important Safety Information			
Parts Detail and Product Specifications			3
Cushion Setup			3-4
Cover Removal and Replacement			5
Cushion Insert Removal and Replacement			
Foam Insert Removal and Replacement			5-6
Troubleshooting			6
Cleaning and Disinfecting			7-8
Storage, Transport, Disposal and Recycling			8
Limited Warranty			
Español	9	Português	81
Français	17	hrvatski	89
Deutsch	25	Čeština	97
Italiano	33	Íslenska	105
Nederlands	41	Polski	113
Dansk	49	Română	121
Svenska	57	Slovenčina	129
Norsk	65	Slovenščina	137
Suomi	73	Türkçe	145

Symbols Glossary

 Manufacturer	 FITS Fits	 WDH Size (width x depth x height)	 FRONT Front of product
 Date of manufacture	 REF Model number	 Weight limit	 BACK Back of product
 Consult instructions for use	 SN Serial number	 MD Medical Device	 BTM Bottom of product
 QTY Quantity		 EC REP Authorized representative in the European Community	

Please Note

Supplier: This manual must be given to the user of this product.
 Operator (Individual or Caregiver): Before using this product, read the instructions and save for future reference.
 Products may be covered by one or more U.S. and foreign patents and trademarks, including DRY FLOATATION®, Hybrid Elite®, Hybrid Elite SR™, ROHO®, Sensor Ready®, shape fitting technology®, Smart Check® and Yellow Rope®. Zytel® is a registered trademark of DuPont.
 Product information is changed as needed; current product information is available at permobilus.com.
 ROHO, Inc. is part of Permobil Seating and Positioning.
 © 2009, 2020 Permobil
 Rev. 2020-04-15

The following products are featured in this operation manual:
 Hybrid Elite Dual Compartment Cushion
 Hybrid Elite SR

Intended Use

The ROHO Hybrid Elite Cushion (Hybrid Elite Cushion) is an adjustable, air-filled, wheelchair support surface combining a ROHO DRY FLOATATION cellular design cushion insert with a contoured foam base and is intended to conform to an individual's seated shape in the pelvic loading area to provide skin/soft tissue protection, positioning, and an environment to facilitate wound healing. There is a weight limit of 500 lbs. (225 kg), and the cushion must be properly sized to the individual. The Hybrid Elite Cushion must be used with the supplied cover. Cushion models include:

- ROHO Hybrid Elite SR Cushion (Hybrid Elite SR)
- ROHO Hybrid Elite Dual Compartment Cushion (Hybrid Elite Dual Compartment Cushion)

Contraindication: Based on clinical, scientific or engineering evidence, the Hybrid Elite Dual Compartment Cushion may not be suitable for individuals who require accommodation of pelvic asymmetry greater than 2 in. (5cm). The Hybrid Elite SR Cushion is not suitable for individuals who require accommodation of pelvic asymmetry.

Compatibility: Hybrid Elite Cushions are NOT COMPATIBLE with the ROHO Contour Base.

A clinician who is experienced in seating and positioning should determine whether the product is appropriate for the individual's particular seating needs. The product's intended use statement alone should not be relied upon to make that determination.

The medical products featured in this document are intended to be only one part of an overall care regimen that includes all seating and mobility equipment and therapeutic interventions. A clinician who is experienced in seating and positioning should identify this care regimen after evaluating the individual's physical needs and overall medical condition.

A clinician should also evaluate for visual, reading, and cognitive impairment to determine the need for caregiver assistance or the need for other assistive technology, such as large-print instructions, to ensure proper product use.

Important Safety Information

Warnings:

- Skin/soft tissue breakdown can occur due to a number of factors, which vary by individual. Check skin frequently, at least once a day. Redness, bruising, or darker areas (when compared to normal skin) may indicate superficial or deep tissue injury and should be addressed. If there is any discoloration to skin/soft tissue, STOP USE immediately. If the discoloration does not disappear within 30 minutes after disuse, immediately consult a healthcare professional.
- DO NOT use the cushion on top of, or in conjunction with, any other products or materials, except otherwise indicated in this manual. Doing so may cause the individual to become unstable and vulnerable to falling. If your cushion feels unstable, see "Troubleshooting" or consult with a clinician for proper use.
- DO NOT place obstructions between the individual and the cushion. Except for the compatible covers and accessories listed in "Product Specifications" this manual, any items placed between the individual and the cushion, whether inside or outside the cover, will make the cushion less effective.
- The cushion and the cover MUST be compatible sizes and MUST be used as directed in this manual. IF THEY ARE NOT: 1) the benefits provided by the cushion may be reduced or eliminated, increasing the risk to skin or to other soft tissue, and 2) the individual may become unstable and vulnerable to falling.
- DO NOT expose the product to high heat, open flames or hot ashes. Testing or certification claims, including for flammability, may no longer apply to this device when it is combined with other products or materials. Check testing and certification claims for all products used in combination with this device.
- The cushion surface will acclimate to the surrounding temperature. Take appropriate precautions, especially when the cushion will contact unprotected skin.

Warnings:

- If you are unable to perform any tasks described in this manual, seek assistance by contacting your clinician, equipment provider, distributor, or Customer Support.
- DO NOT use the components of this product separately. The cushion must be used with all components assembled.
- DO NOT use the cushion as a water flotation device (e.g., a Life Preserver). It will NOT support you in water.

Cautions:

- Keep the cushion away from sharp objects.
- Changes in altitude may require adjustment to the cushion.
- If the cushion has been in temperatures lower than 32° F (0° C) and exhibits unusual stiffness, allow the cushion to warm to room temperature.
- DO NOT use a pump, cover or repair kit other than those compatible with your product
- DO NOT modify the cushion or any components. Doing so may result in damage to the product and will void the warranty.
- DO NOT allow your cushion to come into contact with oil-based lotions or lanolin. They may degrade the material.
- Prolonged exposure to ozone may degrade materials used in the cushion, affect the performance of the cushion, and void the product warranty.
- Periodically, check components for damage and replace if necessary. See "Troubleshooting".
- DO NOT use the inflation valve or quick disconnect as a handle for carrying or pulling the cushion. Carry the cushion by the cushion base or by the carrying handle on the cover.

Parts Detail and Product Specifications

Package contents: Hybrid Elite cushion, cover*; repair kit, hand pump, operation manual, limited warranty supplement, product registration card.

* Instructions are provided with the cover.

To order replacement parts, contact your product provider or Customer Support.

The serial number label is located on the bottom of the cushion.

Materials: Cushion Insert: neoprene rubber, made without latex; Inflation Valve: Zytel; Quick Disconnect: acetal; Contoured Foam Base: polyethylene closed cell foam; Foam Insert: polyurethane foam

Cover materials - Refer to the care label on the cover.

Cushion expected lifetime: 5 years

Covers and Accessories Compatible with Hybrid Elite Cushions

Note: Not available in all markets. Contact Customer Support.

ROHO Heavy Duty Cushion Cover


ROHO Hybrid Elite Cushion Cover

ROHO Planar Solid Seat Insert

ROHO Repair Kit

ROHO Hand Pump


Smart Check Cushion Monitoring Device (Hybrid Elite SR)

A: Dual Inflation Valve (Hybrid Elite Dual Compartment Cushion)

B: Quick Disconnect (Hybrid Elite SR)

WIDTH (W) - in. (cm) ⇔	13¾ (35)	14¾ (37.5)	15¾ (40)	16¾ (43)	17¾ (45)	19¾ (50)	21¾ (55.5)	23¾ (60.5)
DEPTH (D)** - in. (cm) ⇓	Contoured Foam Base width x depth: 1RHExxxxC-SR & 2RHExxxxC*							
14¾ (37.5)	1414							
15¾ (40)	1415	1515	1615					
16¾ (43)	1416	1516	1616	1716	1816			
17¾ (45)		1517	1617	1717	1817	2017		
18¾ (48)			1618	1718	1818	2018	2218	2418
20¾ (53)			1620	1720	1820	2020	2220	2420
22¾ (58)							2222	

Note: All measurements are approximate. Width and depth measurements are taken from the top of the foam base and are rounded to the nearest ¼ in. or 0.5 cm. Cushion item numbers and cover tags indicate the wheelchair size the cushion is designed to fit.

* Indicates model; for example, 1RHE1616C-SR or 2RHE1616C

Weight: 3¼ lbs. / 1.4 kg (based on 1RHE1817C-SR) Height: 4 in. / 10.0 cm

** Depths shown include ¾ in. (1.5 cm) of foam base that fits between the wheelchair back posts.

Cushion expected lifetime: 5 years

IF YOU OWN A SMART CHECK: The Hybrid Elite SR is equipped with Sensor Ready Technology and can be used with a Smart Check, which is a cushion feedback system. Smart Check can be used to find the individual's inflation level and help check their cushion. For more information, contact Customer Support.

Note: If you own a Smart Check, STOP HERE and refer to the setup instructions provided with the Smart Check.

Cushion Setup

⚠ Warnings:

- DO NOT use a product that is underinflated or overinflated, because 1) the product benefits will be reduced or eliminated, resulting in an increased risk to skin and other soft tissue, and 2) the individual may become unstable and vulnerable to falling. Carefully follow the instructions for inflation, placement and hand check. If the product does not appear to be holding air, or if you are not able to inflate or deflate the product, see "Troubleshooting". Immediately contact your equipment provider, distributor, or Customer Support if the problem persists.
- When possible, use hook and loop fasteners to further secure the cushion to the seating surface. Hook fasteners are available upon request.
- After setting up the cushion the first time, perform a hand check frequently, at least once a day. When the cushion is used by a different individual or in a different wheelchair, repeat cushion setup. Follow the hand check instructions in this manual.

Cushion Setup is continued on the next page.

Cushion Setup, continued


! **IMPORTANT!** For dual compartment cushions, repeat all steps for both valves.

Inflate, cover, and transfer.


1 Inflate: Open the cover flap. For the Heavy Duty Cover, open the zipper.

- Slide the hand pump nozzle over the inflation valve. Open the inflation valve. Inflate the cushion until all of the air cells feel firm.
- Close the inflation valve until it comes to a stop. Remove the hand pump. Make sure the inflation valve is closed.

Note: To open, turn the inflation valve at least one full rotation. The inflation valve will spin freely when open.

Open

Close

2 Cover: Make sure the cover is correctly installed. Refer to the instructions provided with the cover and the "Cover Removal and Replacement" section in this manual.

3 Transfer:

- Place the cushion on the wheelchair with the air cells facing up and the non-skid surface of the cover facing down. The symbol on the cover showing the seated figure may be used to help orient cushion. (Fig.1)
- The individual should transfer to the wheelchair and sit in their usual sitting position. Confirm that the bony prominences are supported by the air cells. The inflation valve should be in the front, on the left side of the seated individual. (Fig. 2)


Fig. 1

Fig. 2

Perform a hand check.

4 To perform a hand check, use the width of your fingers to check the amount of air in the cushion. When the cushion has the correct amount of air*, the individual is both immersed in the cushion and supported by air:

- Slide your hand between the air cells and the individual, and feel for the lowest bony prominences. *Suggestion: To help locate the bony prominences while performing a hand check, lift and lower the leg.*
- Refer to the Hand Check Guide below to adjust the cushion.


* The recommended distance between the lowest bony prominence and the hard surface is ½ - 1 in. (1.5 - 2.5 cm).

HAND CHECK GUIDE

! **IMPORTANT!** For dual compartment cushions, check both compartments.


UNDERINFLATED!

Not able to move fingers.
Not enough air in the cushion.
Risk: Not supported by air.


OVERINFLATED!

Can move fingers too freely.
Too much air in the cushion.
Risk: Not immersed in the cushion.

CORRECT AMOUNT OF AIR!

Able to slightly move fingers.
Immersed in the cushion and supported by air.

If the cushion is UNDERINFLATED:

Add some air. Close the inflation valve. Perform a hand check again. Repeat until the cushion has the correct amount of air.

If the cushion is OVERINFLATED:

Remove some air. Close the inflation valve. Perform a hand check again. Repeat until the cushion has the correct amount of air.

If your cushion has the CORRECT AMOUNT OF AIR, it is ready to use.

! **Warning:** Make sure there is at least ½ in. (1.5cm) of air supporting the individual's lowest bony prominences.

Cover Removal and Replacement

IMPORTANT: To remove or replace the cover, refer to the instructions provided with the cover and follow these warnings and the cover flap instructions.


Warnings:

- Before replacing the cover, make sure the foam insert and the cushion insert are correctly positioned and securely attached. See "Foam Insert Replacement" and "Cushion Insert Replacement".
- If hook fastener strips were removed prior to washing, replace the strips.

Cover Flap Instructions (Hybrid Elite Cover only; does not apply to the Heavy Duty Cover)

Removing the Hybrid Elite Cover: With the cover unzipped, reach inside the cover and pull the inflation valve through the opening so that it is inside the cover before removing the cover from the cushion. Repeat for both sides.


Installing the Hybrid Elite Cover: With the cover unzipped, reach inside the cover and push the inflation valve through the opening beneath the cover flap, and then close the cover flap. Repeat for both sides.

Cushion Insert Removal and Replacement


Cushion Insert Removal

- 1 Remove the cover from the cushion. Refer to the instructions provided with the cover and to the "Cover Removal and Replacement" section in this manual.
- 2 Slide the tubing out of the loop on the foam insert.
- 3 Remove the tubing from the groove in the foam base and gently pull through the hole in the foam base corner. Repeat for other side.
- 4 Remove the cushion insert from the contoured foam base.

Cushion Insert Replacement


- 1 Place the cushion insert on the foam base, lining up the hook and loop fastener strips. Make sure the cushion insert is centered. Press down to securely attach the cushion insert to the foam base.
- 2 Insert the tubing through the hole in the foam base corner. Repeat for other side.
- 3 Slide the tubing through the loop on the foam insert and press the tubing into the groove in the foam base. Repeat for other side.
- 4 Confirm correct positioning of the foam insert. See "Foam Insert Replacement".
- 5 Replace the cover. Refer to the instructions provided with the cover and to the "Cover Removal and Replacement" section in this manual.

Foam Insert Removal and Replacement

Foam Insert Removal

- 1 Remove the cover from the cushion. Refer to the instructions provided with the cover and to the "Cover Removal and Replacement" section in this manual.
- 2 Slide the tubing out of the loop on the foam insert. Repeat for other side.
- 3 Flip the cushion over so that the cells face down and the bottom of the contoured foam base faces up.
- 4 Pull the foam insert flap away from the hook fastener pieces on the bottom of the contoured foam base. Repeat for the other side.
- 5 Flip the cushion over so that the cells face up and the bottom of the contoured foam base faces down. Release the foam insert from the contoured foam base by pulling away from the hook fastener pieces.

See the next page for Foam Insert Replacement.




Foam Insert Removal and Replacement, continued

Foam Insert Replacement

Warning: The foam insert must be correctly positioned on the contoured foam base.

Note: If the cushion insert was removed from the contoured foam base, replace the cushion insert before replacing the foam insert. See "Cushion Insert Replacement".

- 1 Position the foam insert on the contoured foam base.
 - The foam insert should be next to the air cells (A). Make sure the foam insert is not on top of the air cells.
 - The foam insert should line up with the front edge of the foam base (B). Make sure the foam insert does not extend past the front edge of the foam base.
 - The foam insert should be centered on the foam base.
- 2 Once correctly positioned, press down on the foam insert to securely attach it to the contoured foam base.
- 3 Slide the tubing through the loop on the foam insert. As needed, push the tubing back into the groove. Repeat for other side.
- 4 Flip the cushion over so that the cells face down and the bottom of the contoured foam base faces up.
- 5 Pull and attach the foam insert flaps to the hook fastener pieces on the bottom of the contoured foam base.
- 6 Flip the cushion over so that the cells face up and the bottom of the contoured foam base faces down.
- 7 Replace the cover. Refer to the instructions provided with the cover and to the "Cover Removal and Replacement" section in this manual.

Troubleshooting

For additional assistance, contact your equipment provider, distributor, or Customer Support.

The cushion is not holding air.	Inflate the cushion. Inspect for damage to inflation valve(s), hose(s) and the quick disconnect (Hybrid Elite SR). Confirm that the inflation valve(s) are completely closed. Look for holes in the cushion. If very small holes or no holes are visible, follow the instructions in the Repair Kit provided with your product. For damage to the inflation valve(s), hose(s), or quick disconnect, or for large holes or leaks in the cushion, see the Limited Warranty supplement, or contact Customer Support.
Not able to inflate or deflate the cushion.	Confirm that the inflation valve(s) are open. Inspect the inflation valve(s) and hose(s) for damage. For the Hybrid Elite SR: also check the quick disconnect for damage.
The cushion is uncomfortable or unstable.	<ul style="list-style-type: none"> - Make sure that the air cells are not underinflated or overinflated (see "Cushion Setup"). - Make sure that the foam insert and the cushion insert are correctly placed on the contoured foam base. - Make sure that the foam insert does not overlap the air cells. - The air cells and the top of the cover should face up. - The cushion should be centered beneath the individual. - The cushion should be properly sized to the individual and the wheelchair. Allow the individual at least one hour to adjust to the cushion.
The cushion slides on the wheelchair.	Make sure that the non-skid bottom of the cushion cover faces away from the individual. Use directional labels on the cushion and on the cover to orient the cushion. When possible, use hook and loop fasteners to further secure the cushion to the wheelchair. Hook fasteners are available upon request. Refer to the instructions provided with the cover.
A component is damaged.	Contact Customer Support for a replacement.

Cleaning and Disinfecting

Cover Care

For cover care instructions, refer to the manual provided with the cover.

Cushion Care

Warnings:

- DO NOT machine wash or machine dry the contoured foam base or the foam insert. DO NOT machine dry the cushion insert.
- DO NOT allow water or cleaning solution to enter the cushion. Make sure the inflation valve(s) are closed. Hybrid Elite SR: Make sure the Smart Check is detached from the cushion. Water in a cushion with Sensor Ready (SR) technology can result in damage to the Smart Check and cause a false reading when using Smart Check.
- Cleaning and disinfecting are separate processes. Cleaning must precede disinfection. If the product becomes soiled, or between use by different individuals: clean, disinfect, and check the product for proper functioning.
- Clean the product regularly and when it becomes soiled.
- Follow all manufacturer bleach container and germicidal disinfectant safety guidelines.
- Note about sterilization: High temperatures accelerate aging and will damage the product assembly. The products featured in this manual are not packaged in a sterile condition, nor are they intended or required to be sterilized prior to use. If institution protocol demands sterilization: First, follow the cleaning and disinfecting instructions; then, open the inflation valve(s) and use the lowest possible sterilizing temperature, but not to exceed 70°C (158°F), for the shortest time possible. DO NOT autoclave.

Cautions:

- Use of the following may cause damage to the cushion: abrasives (steel wool, scouring pads); caustic, automatic dishwasher detergents; cleaning products that contain petroleum or organic solvents, including acetone, toluene, Methyl Ethyl Ketone (MEK), naphtha, dry cleaning fluid, adhesive remover; oil-based lotions, lanolin, ozone gas; ultraviolet light cleaning methods.
- DO NOT expose to direct sunlight.

Recommendation: When washing by hand, use room-temperature water.

Note: Deflated air cells may be easier to thoroughly clean. To deflate the air cells: with the inflation valve(s) open, roll up and/or compress the cushion, and then completely close the inflation valve(s).

Hand Wash and Disinfect the Cushion Insert (Preferred Method)

Carefully remove the cushion insert from the contoured foam base. Close the inflation valve(s). Hand wash, using mild liquid hand soap, hand dishwashing liquid, laundry detergent, or a multipurpose detergent (Follow product label instructions.) Gently scrub all cushion surfaces, using a soft plastic bristle brush, a sponge or a cloth. Rinse with clean water. Air dry.

To disinfect, repeat cleaning instructions, using 1 part household liquid bleach per 9 parts water. Keep the product wet with the bleach solution for 10 minutes. Rinse thoroughly with clean water. Air dry.

 **Caution:** Rinse thoroughly. Detergent residue may cause air cells to stick together.

Clean and Disinfect the Contoured Foam Base and Foam Insert

To clean and disinfect the Contoured Foam Base: Remove the cover, foam insert and cushion insert. To clean, wipe gently with a cloth dampened with water and a household detergent, then rinse with cool, clean water.

To disinfect, follow the cleaning instructions, but use disinfectant instead of household detergent. After cleaning or disinfecting, allow the foam to air dry thoroughly before putting the cover on.

To clean and disinfect the Foam Insert: To clean, moisten a cloth, then wring to remove excess water. Wipe gently with the cloth and a household detergent, then wipe again with a clean, slightly damp cloth.

To disinfect, follow the cleaning instructions, but use disinfectant instead of household detergent. After cleaning or disinfecting, allow the foam to air dry thoroughly before putting the cover on.

Cleaning and Disinfecting, continued

Machine Wash and Disinfect the Cushion

⚠ IMPORTANT: Machine washing can cause wear and damage, including cuts and tears to the cushion and the detachment of required product information. Washing and disinfecting the cushion by hand and air drying is strongly recommended.

⚠ Warnings:

- DO NOT machine wash or machine dry the contoured foam base or the foam insert. DO NOT machine dry the cushion insert.
- DO NOT allow water or cleaning solution to enter the cushion. Make sure the inflation valve(s) are closed. Hybrid Elite SR: Make sure the Smart Check is detached from the cushion. Water in a cushion with Sensor Ready (SR) technology can result in damage to the Smart Check and cause a false reading when using Smart Check.
- High water temperature may cause cushion discoloration and label detachment.
- DO NOT expose cushions to temperatures above 70°C (158°F), which may damage cushion components.
- For machine washing or disinfecting, make sure the washing machine capacity is sufficient. A washing machine that isn't large enough may not remove all soil or contaminants from a cushion.
- Using a damaged cushion will reduce or eliminate the benefits of the cushion and increase risk to skin and other soft tissue.

Prepare the Cushion for Machine Cleaning:

- 1) Remove the cover and close the inflation valve(s). Make sure the Smart Check is detached from the cushion (Hybrid Elite SR). Examine the cushion for damage, and patch any holes (refer to Repair Kit instructions provided with the product).
- 2) Remove stains or odors as needed: Place the cushion in warm water and a disinfectant cleaning agent. Gently scrub, using a soft plastic bristle brush, and rinse with clean water.
- 3) Deflate the cushion: Open the inflation valve and roll up the cushion, starting at the corner opposite the valve. After you have removed as much air as possible, close the inflation valve.
- 4) Inspect for any sharp edges in the washing machine that may damage the cushion.

Machine Wash the Cushion:

Place the cushion in a laundry bag and then into the washing machine. When folding the cushion, make sure the air cells face outward. Use warm water (85° - 105°F / 30° - 40°C) and laundry detergent or machine- and cushion-safe cleaner/biocide, placed in the detergent dispenser. Gentle cycle: Use a machine setting that does not include a high-speed spin cycle.

Machine Disinfect the Cushion: Follow the machine washing instructions and use hot water 60°C (140°F). Use biocide (refer to the manufacturer instructions). Bleach with care: Only use a dilution of 1 part household liquid bleach per 9 parts water. Rinse with clean water.

⚠ Caution: Rinse thoroughly. Detergent residue may cause air cells to stick together.

Drying the Cushion: Remove the cushion from the laundry bag and air dry. For faster drying, use a hand-held hair dryer, hold no closer than 6 in. (15 cm) from the cushion, and use a sweeping motion.

⚠ Warning: Overheating the neoprene rubber can cause permanent product damage.

Care Symbols

Hand wash.

Line dry in the shade.

Bleach (1 part bleach: 9 parts water).

Machine wash, warm, gentle.

Machine wash hot, normal, at the temperature shown.

Use a handheld dryer.

Wipe with a damp cloth.

Hand wash cold.

Dry flat in the shade.

Storage, Transport, Disposal and Recycling

Storage and Transport: Clean and disinfect the product before storing. Open the inflation valve(s). Deflate the product. Store the product in a container that protects the product from moisture, contaminants and damage. After storing, follow instructions to prepare the product for use. After transport, inspect the product to ensure it was not damaged (see Troubleshooting).

Disposal: The components of the products in this manual are associated with no known environmental hazards when properly used and when disposed of in accordance with all local/regional regulations. At the end of its useful life, treat as healthcare waste and dispose of according to local healthcare-waste guidelines. Incinerate with other healthcare waste. Incineration must be performed by a qualified, licensed waste management facility.

Recycling: Contact a local recycling agency to determine recycling options for the product.

Limited Warranty

Limited Warranty Term: 24 months from the date the product was originally purchased. The warranty does not apply to punctures, tears, burns, or misuse. See also the instructions provided with the cover and the Limited Warranty supplement provided with your product, or contact Customer Support.

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ROHO, Inc. pertenece Permobil Seating and Positioning.

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